



Arm Sling with Padded Shoulder

APPLICATION INSTRUCTIONS:

1. Unfold the sling.
2. Insert hand and lower arm into sling in the appropriate position.
3. Disconnect the Velcro.
4. Bring the end of the strap through the D-ring, adjust to fit and re-apply the Velcro.
5. Adjust the neck pad to comfort.
6. (Option) The thumb rest loop on the inside of the sling is optional.

Washing Instructions:

Hand wash in cold water using mild soap. Do not bleach. Line dry only. Avoid direct heat or sun.

NOTE: Do not over tighten, adjust as necessary.

CAUTION: If swelling increases and pain, numbness, skin irritation occurs, discontinue use. Contact medical provider immediately.

This product is only a support device, not to prevent injuries or prevent reoccurrence of injuries.



OvationMedical[®] 

5155 Clareton Dr.
Agoura Hills, CA. 91301

USA

800 403-6466

www.ovationmed.com

MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany