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






ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS: The Pneumatic Short Walking Brace is intended for use as a support for lower leg, ankle, and toe region, as a fixation of foot and ankle at 90 degrees that allows ambulation. These products are indicated for use by individuals suffering from any one of the following, or other similar conditions; lower leg fractures, fractures of the metatarsals, ankle sprains, ligament injuries, post fractures after initial casting, post-operative procedures, among other musculoskeletal ailments typically seen in an orthopedic clinic or hospital. The Pneumatic Short Walking Brace is specifically indicated for soft tissue injuries, Grade 2 and 3 sprains-stable fracture, post-operative use, trauma and rehab.

The Pneumatic Short Walking Brace is intended for short term single patient use and not intended to undergo regular service cycles.

This product was designed to withstand normal use conditions for orthopedic support, post injury rehab or prevention of further injury. Like all products that are worn by individuals, the extended life of such products is determined by the level of activity and type of activity of the individual.

CONTRAINDICATIONS: None known. In every case, the final decision whether to use this Ovation Medical® device rests with the attending physician/therapist. He/she should assess prior to use of the device, if it is suitable for a given application, taking into account the condition of the patient.

WARNINGS AND PRECAUTIONS: Do not over-inflate aircells. Improper aircell inflation may cause significant skin irritation, especially in people with diabetic neuropathy (a diabetic nerve condition). Check your skin frequently for "hot spots" and skin irritation. If you see any skin changes or feelings of discomfort, reduce the aircell compression. Do not use this device if you cannot feel sensations whether due to post-op anesthesia, nerve, skin or other conditions.

This product is not intended to be used directly on open wounds or broken skin. If pain and/or numbness occurs, contact your medical professional immediately. Ask prescribing medical professional as to what opposing footwear to use.

Use caution when walking on wet or slippery surfaces to avoid injury. At high altitudes, the aircells will expand beyond optimal level. Adjust aircell pressure by tightening or loosening straps.

Always consult with your physician or therapist before making changes to the brace. This product is to be used under the supervision of a medical professional. This device is not intended for unsupervised public use.

REPORTING PROBLEMS: The user should report any suspected serious incident related to the use of this Ovation Medical® device by informing the manufacturer and the competent authority of the state or the country in which the incident occurred.

APPLICATION INSTRUCTIONS:

- Unfasten Velcro straps and remove inner soft good liner from boot.
- Open the soft good liner and place foot inside with the heel fitting snugly into the back portion of the liner. Fasten and wrap foot flaps then the leg portion of the liner. Ensure liner is snug but not constrictive. Make sure pump is facing outward.
- Gently spread the side struts and place the lower leg into the brace.
- The ankle should be centered at the struts with the foot and ankle at 90°.
- Fasten the Velcro straps beginning at the toe then moving up the leg. Do not over tighten. As swelling pattern changes throughout the day*, it may be necessary to adjust the straps for comfort. *Due to activities, rest, etc.
- To inflate air liner, press and release center air bulb until liner provides a comfortable cushion and support. Do NOT over inflate. To deflate or adjust compression, press down on smaller circular side button.
- To remove, simply unfasten the straps without removing them completely from the assembled walker. Unfasten the soft good liner and gently bring the foot and ankle forward and out of the walker. To reapply, place the foot and ankle into the liner and reattach the Velcro closures (As per Fig. 5).

CARE AND MAINTENANCE: Hand wash in cold water using mild soap, do not bleach, line dry only, avoid direct heat or sun. Store in a cool dry place.

DEVICE DISPOSAL: Ovation Medical® products should be properly disposed of by a qualified recycler when the post-operative use period is completed or when the device is not longer required by the patient.

SYMBOLGY:

- Consult instructions for use and accompanying documents.
- Manufactured by. Note: can be used in combination with date of manufacture.
- Date of manufacture accompanied by a date.
- Single patient multiple use.
- Device sold in non-sterile condition.
- Not made with Natural Rubber Latex
- Peel here.
- Lot number.
- Reference number or product number.
- Quantity contained in package.
- UDI Unique Device Identifier Number.
- Product intended to be used as a medical device.
- Authorized representative (REP) in the European Community (EC).
- Conform to the European Medical Device Regulation (EU) 2017/745 for class I medical devices

Patents: 9492305, 10039664, 10085871, D759953. Additional patents applied for.

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