



Shoulder Abduction Sling

APPLICATION INSTRUCTIONS:

1. Apply the Sling portion first.
2. Unfasten the Velcro at the top of the sling
3. Insert hand and lower arm into sling in the appropriate position.
4. Bring the end of the strap through the D-ring, adjust to fit and re-apply the Velcro.
5. Adjust the neck pad to comfort.
6. (Option) The thumb rest loop on the inside of the sling is optional.
7. Apply pillow in comfortable position behind the arm sling and with the curved side against the body.
8. Attach pillow Velcro to arm sling Velcro.
9. Adjust waist strap as necessary and thread the end through the d-ring. The strap may be cut to modify positioning of the Velcro.

Washing Instructions:

Hand wash in cold water using mild soap. Do not bleach. Line dry only. Avoid direct heat or sun.

Caution: In the event of discomfort or numbness of any kind, consult your caregiver immediately.

NOTE: Do not over tighten, adjust as necessary.

CAUTION: If swelling increases and pain, numbness, skin irritation occurs, discontinue use. Contact medical provider immediately.

This product is only a support device, not to prevent injuries or prevent reoccurrence of injuries.



OvationMedical® 

5155 Clareton Dr.
Agoura Hills, CA. 91301
USA
800 403-6466
www.ovationmed.com

MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany